

C-Chip

CE

The World's First True Disposable Hemocytometer



Precise volume control
Easy to use &

NanoEnTek

e-mail
website

sales@nanoentek.com
www.nanoentek.com

NanoEnTek Inc.

12F, Ace High-end Tower, 235-2, Guro3-dong,
Guro-gu, Seoul, 152-740, Korea
Tel : +82-2-6220-7913 Fax : +82-2-6220-7721

NanoEnTek USA Inc.

5627 Stoneridge Drive Suite 304, Pleasanton,
CA 94588, USA
Tel : +1-925-225-0108, +1-888-988-0108(Toll free)
Fax : +1-925-225-0109

EC REP

MT Promedt Consulting GmbH
Altenhofstrasse 80,
66386 St. Ingbert, Germany

C- Chip is

a Disposable Plastic Hemocytometer

- Does not need coverslips
- Eliminates the need of washing and reuse
- Reduces exposures to the infectious samples and contamination
- Guaranteed repeatability and reliability
- Two test
- Quartz grade optical plastic
- Sturdy and strong






Applications

Blood analysis (Hematology) | Blood cell counting

Cell Culture | Cell concentration measurement / Cell viability test

Microbiology | Bacteria & Fungal spore counting

IVF, IUI | Sperm counting



Models

No.	Grid	Depth
DHC-N01	Neubauer Improved	100 μm
DHC-F01	Fuchs Rosenthal	200 μm
DHC-B01	Burker	100 μm
DHC-B02	Burker Turk	100 μm
DHC-M01	Malassez	200 μm

Specifications

Dimension : 25 (W) x 75 (L) x Thickness 1.6 mm

Chamber depth : 100, 200, 10 μm

Chamber volume : 10 μL

NanoEnTek

e-mail sales@nanoentek.com
website www.nanoentek.com

NanoEnTek Inc.

12F, Ace High-end Tower, 235-2, Guro3-dong,
Guro-gu, Seoul, 152-740, Korea
Tel : +82-2-6220-7913 Fax : +82-2-6220-7721

NanoEnTek USA Inc.

5627 Stoneridge Drive Suite 304, Pleasanton,
CA 94588, USA
Tel : +1-925-225-0108, +1-888-988-0108(Toll free)
Fax : +1-925-225-0109

EC REP

MT Promed Consulting GmbH
Altenhofstrasse 80,
66386 St. Ingbert, Germany